

# Using Micro Computed Tomography to Better Characterize the Apollo Sample Suite: a Retroactive PET Style Analysis.



Ryan Zeigler: NASA – Johnson Space Center

Jeremy Kent: GeoControl Systems, Jacobs JETS Contract

Charis Krysher: HX5, Jacobs JETS Contract

Andrea Mosie: GeoControl Systems, Jacobs JETS Contract

Linda Watts: Barrios Technology, Jacobs JETS Contract



# Motivation

- The overarching goal of this presentation is to show that the micro-CT technique is a useful tool in the curation of Apollo samples.
  - Most of what I will show here is applicable to other astromaterials collections at JSC (esp. meteorites)
- Ultimately we used this as the justification for adding micro-CT as part of the core capabilities of JSC curation.
- Part of our larger plan to modernize the curation process for Apollo (and all other) samples, using non- (or minimally) invasive techniques to improve the curation process.
  - Would include other techniques in the future, e.g., micro-XRF, scanning Raman, FTIR, 3-D scanning and imaging, etc.

# More Motivation

- The Apollo sample collection is a “mature” collection
  - Most of the large samples have been well studied, though new advances are still being made with them.
- PIs continue to clamor for new samples to study
  - Plutonic samples to nail down the timing of the formation and evolution of the Moon.
  - Evolved lithologies to constrain how the lunar crust evolved.
  - Better understand the volatile budget of the Moon
- Terrestrial weathering in lunar meteorites limits their usefulness for some things.
- The key to giving PIs what they want can be found in clasts in polymict Apollo breccias
  - We need to characterize these while preserving their pristinity

# Apollo Sample Suite

- From 1969-1972 the Apollo astronauts collected 382 kg of rocks, soils, and cores from 6 locations on the Moon.
- They are the only astro-material samples collected with geologic context, which increases their usefulness.
- Studies of the large igneous rocks from the Moon shaped (and continue to shape) our understanding of lunar formation and evolution.


76535  
Troctolite  
Mg-Suite



62275 - FAN



15556  
Low-Ti  
Basalt






# Apollo Sample Suite

- The majority of Apollo rock samples returned are in fact polymict breccias.
  - If we exclude basalts it is an overwhelming majority.
- These contain thousands of small (mm to cm) sized clasts.
  - Modern analytical techniques have reduced the mass needed for analyses.
- Unambiguously identifying the lithology of a clast using optical microscopy is challenging.
- Many polymict breccias are glass coated, further obfuscating classification of clasts.



# Micro-CT – The Basics

- High energy (up to  $\sim 450$  kV, 100-200 kV typical for hard rock scans) x-ray beam is projected through a sample.
  - Uses a micro-focus x-ray source (typically  $1\text{-}5\text{ }\mu\text{m}$ )
- The attenuation of the X-ray cone is recorded on a detector in 2-D slices
  - Many 2-D slices are acquired as the sample is rotated by fractions of a degree ( $>1000$ )
  - The resolution of the scan is related to the size of the sample and the detector.
- Computer software is used to reconstruct these 2-d slices into 3-d volumes.






# Our Analyses

- We've scanned 10 samples at the Natural History Museum, London
  - 10057,19 – Mare Basalt, 132 g
  - 14305,483 – Regolith Breccias, 155 g
  - 14321,1404 – Regolith Breccia, 346 g
  - 15205,0 – Regolith Breccia, 136 g
  - 15405,0 – Impact Melt Breccia, 201 g
  - 60639,0 – Regolith Breccia, 153 g
  - 15555,62 + ,1039 – Mare Basalt, 33/66 g
  - 15556,111 +, 215 – Mare Basalt, 44/354 g
  - 60015,180 – Ferroan Anorthosite, 122 g
  - 70017,72 – Mare Basalt, 127 g
- Scans were made at 205-220 kV, 135-160  $\mu$ A, and 20-40  $\mu$ m/voxel

Special thank you to Univ. of Texas CT facility as well.



Nikon X-Tek HMXST 225 system at NHM - London




14321 sample triple bagged in Teflon and immobilized in a cardboard tube

# Scan of 60639


- Apollo 16 regolith breccia 60639 is a fairly generic Apollo 16 regolith breccia, with one important caveat:
  - It is the only Apollo 16 sample containing a sizable basalt clast, which is the only mid-TiO<sub>2</sub> basalt on the Moon (6 wt% TiO<sub>2</sub>).









10.00 mm



Subvolume scan –  
20  $\mu\text{m}$  resolution

60639,0  
Apollo 16 polymict breccia

View of the clasts within the sample.




Basalt clast  
segmented out

Micro-CT scan at  
215 kV, 160  $\mu\text{A}$ , 3100 scans

# Scan of 14321


- Apollo 14 regolith breccia 14321 was the largest Apollo 14 sample (> 9 kg) containing important clast types
  - These include troctolites, high-Al and high-K basalt clasts, as well as a variety of evolved clasts such as granites and alkali anorthosites.






# Scan of 14321

- Apollo 14 regolith breccia 14321 was the largest Apollo 14 sample (> 9 kg) containing important clast types
  - These include troctolites, high-Al and high-K basalt clasts, as well as a variety of evolved clasts such as granites and alkali anorthosites.



Relative brightness is caused by the density of the phase and average Z.





MuCalc tool - UTCT




# Scan of 14321

- Apollo 14 regolith breccia 14321 was the largest Apollo 14 sample ( $> 9$  kg) containing important clast types
  - These include troctolites, high-Al and high-K basalt clasts, as well as a variety of evolved clasts such as granites and alkali anorthosites.



Clast =  $0.637 \text{ cm}^3$

Clast =  $\sim 1.9 \text{ g}$



# Other lunar uses for micro-CT...


- In addition to clast identification, micro-CT scans would be useful to address a variety of other lunar sample questions.
  - Identification of mineral grains in igneous lunar samples.
  - Identification of zircons in lunar samples.
  - Measurement of macroscopic porosity in lunar samples



# Other lunar uses for micro-CT...

- In addition to clast identification, micro-CT scans would be useful to address a variety of other lunar sample questions.
  - Identification of mineral grains in igneous lunar samples.
  - Identification of zircons in lunar samples.
  - Measurement of macroscopic porosity in lunar samples


15556,215 – 12 cm tall, 354 g






# Non lunar uses for micro-CT...

- Every single use I just highlighted for lunar samples applies to the Antarctic meteorite collections as well.
  - e.g., CAIs in chondrites, carbonaceous chondrite clasts in Howardites
  - Currently used by the NHM for all martian meteorite allocations.
  - Meteorite classification?
- Micro-CT will likely be a part of the basic curation plan for future lunar sample return missions; it is a base-level requirement for future Mars sample return; and it could be utilized for Osiris Rex.
- It would also be available to internal and external users for purely research purposes






3D SFM Model




Micro-CT "image"



Precision Photography Phase of Lunar Sample 60639 for SFM 3D Modeling. Photographs were taken in 15 degree increments totaling approximately 240 images.

# What would we do once we made the scans?

- A summary of the scan data would be “published” bi-annually to disseminate the information.
  - Similar to the ANSMET newsletter
- Individual slices and higher level data products (e.g., videos) would also be made available to the PIs for additional analysis.
- The scans would be used to identify the best places to make new band saw cuts, exposing the new clasts for more detailed study.
- These newly slabbed surfaces could be better characterized by micro-XRF






RGB = KAlSi

14305,483

Apollo 14 polymict breccia





10 cm

Scans done on the EDAX Orbis PC  
at Washington University

Micro-XRF scans at 40  
kV, 800  $\mu$ A, 30  $\mu$ m beam

RGB = CaMgFe





# Bottom Line

- Adding micro CT as a core part of the JSC curation capabilities will allow us to better utilize our existing collections – Conservation of Mass
- It will allow us to find “new” samples and provide better samples to scientists that allow them to address a wider range of scientific questions.
- At the same time, it will allow us to become proficient in a technique that is going to be important in many future sample return missions.



# 14321 video





Clast = 0.637 cm<sup>3</sup>


Clast = ~1.9 g



# 14321 video




Clast = 0.637 cm<sup>3</sup>  
Clast = ~1.9 g



# Scans of 15405 and 15205

- Pair of Apollo 15 regolith breccias, both of which have been shown to contain abundant evolved and/or KREEPy clasts, e.g., Granites, monzogabbro, etc.





# Micro-CT for Curation

- While it's application to Apollo samples is new, using it on meteorites is not:
  - Numerous scientific studies, e.g., Kuebler et al 1999; Rubin et al 2001; McCoy et al. 2006; etc.
  - Becoming more common in the curation of meteorites, e.g., Smith et al. 2013; Hyde et al. 2013.
- CT imaging of large samples presents a unique set of challenges not necessarily applicable to previous research applications:
  - Resolution is tied to sample size.
  - Sensitivity is also tied to sample size.

**Pebble volumes: 0.731,  
0.190, 0.109, 0.040,  
0.027cm<sup>3</sup>.**

